

SANITECH SERIES 9000 UTENSIL WASHER/DECONTAMINATOR.**INSTALLATION INSTRUCTIONS****Specifications**

Height	1540 mm
Width	600 mm
Depth (closed)	610 mm
Depth (open)	1200 mm

Service requirements

Hot Water	20mm BSP, 50kPa min WM stop tap. Required temperature of hot water supply 60C.
Electrical	240V, 50Hz 10 Amp.
Drain hose outlet	Minimum 38mm standing waste.

INSTALLATION

Before placing the machine in position, prepare the electrical, water and drainage connections. Final connections can then be made when the machine is pushed into position. Hot water washing machine stop tap should be situated above or next to the machine either to the left or right in an accessible position. The machine is supplied with the water inlet pipes already fitted. They are cross-linked polythene pipes (PEX), fitted with a 3/4" BSP internal threaded connector. The pipe can be bent as required: take care, however, not to kink it.

DRAINAGE CONNECTION

The machine is supplied with the outlet hose connected.
Run the outlet hose into top of standing waste.
Drain outlet hose to be 400mm above bottom of door.
Top of standing waste to be minimum of 1000mm from floor.

LEVELING

Adjust leveling feet of the machine to allow unit to sit square.


Push the Machine into Position

When all connections for power supply, water and drainage are ready, push the machine into position, taking up the slack in the outlet hose and supply pipes at the same time. This prevents kinks from being formed in the hose or pipes. Re adjust leveling feet.

SANITECH SERIES 9000 UTENSIL WASHER/DISINFECTOR

Hot water washing machine stop tap should be situated above or next to the machine either to the left or right in an accessible position.

Run the outlet hose into top of 38mm or 50mm standing waste.



STARTING MACHINE FOR THE FIRST TIME

Once the machine has been connected to the power outlet, water connection and turned on:

‘Utensils’ will appear on the Control Panel.

The Control panel will then read ‘D’ and ‘Utensils’.

‘D’ is displayed on the Control Panel indicating that the detergent has not yet been connected.

Connection of detergent is done by referring to page 7 of the Use and Care Guide.

Once the detergent has been connected and primed into the machine the Control Panel will read ‘Utensils’. The machine is ready for normal use.

Please call if you require further consultation.

SANITECH AUSTRALIA

Tel: 61 3 9809 2701

Mob: 0427 823 574

Email: info@sanitech.com.au

Website: www.sanitech.com.au